

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Abstract

Purge and Trap (P&T) is a concentration technique used for the analysis of Volatile Organic Compounds (VOCs). The major component of any P&T system is the analytical trap. This trap is responsible for trapping the analytes purged from the sample and then releasing them upon heating for further analysis by the GC and GC/MS.

This study will evaluate a new analytical trap that is designed by Teledyne Tekmar. The packing material in this new trap improves the recovery and performance of Volatile Hydrocarbon Fractions found in Gasoline Range Organics (GRO) when prepared using Purge & Trap and Headspace sample concentration techniques used in tandem with GC.

Introduction

Purge and Trap (P&T) concentration is a typical technique used for the analysis of Volatile Organic Compounds (VOCs). The major component of any P&T system is the analytical trap. This trap is responsible for retaining the VOCs during purge and then releasing the analytes upon desorbing. The requirements for an analytical trap are as follows:

1. At low temperatures, it must retain the desired analytes while allowing oxygen and water to pass through unimpeded.
2. It must release the analytes quickly and efficiently upon heating
3. Must not contribute any volatiles of interest to the system.
4. Should have a reasonable price and lifetime.

The latest advancement in analytical traps by Teledyne Tekmar is the #11-VPH-Trap. This trap provides the features of Tekmar's #9 analytical trap while modifying the sorbent bed so a lesser amount of methanol is retained. The #11-VPH-Trap is ideal for running the Wisconsin Gasoline Range Organics (GRO) method (e.g. **methyl-tert-butyl-ether, toluene, benzene, ethylbenzene, xylenes, 1, 2, 4-trimethylbenzene, 1, 3, 5-trimethylbenzene and naphthalene**) and BTEX methods.

To evaluate the new #11-VPH-Trap, the analytical trap examined three Tekmar products the Atomx an automated VOC sample prep system (U-Shape), Velocity XPT Purge and Trap (straight trap) with an AQUATEk 100 autosampler, and a HT3™ headspace analyzer (straight trap).

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

For the Atomx, a calibration curve by summation, Replicate Lab Control Spikes (RLCS), Duplicate Lab Control Spikes (DLCS), Limit of Detection (LOD) and Limit of Quantitation (LOQ) were determined for water, high-level soil and low-level soil samples.¹ For the Velocity XPT a calibration by summation, RLCS, DLCS, LOQ and LOD were evaluated for water samples only.¹ The HT3 headspace analyzer generated individual calibration curves for all ten compounds listed above.

Atomx Automated VOC Sample Prep System

Instrument Parameters

The Atomx, equipped with a #11-VPH-Trap, and an Agilent 6890 GC/FID were utilized for this study. **Tables 1-4** show the GC/FID and P&T conditions for water, soil, and methanol extraction applications.

| GC Parameters | |
|---------------|---|
| GC: | Agilent 6890 |
| Detector | FID |
| Column | Restek RTX-VMS 20m x 0.18mmID x 1um |
| Oven Program: | 40°C for 4 min; 16°C/min to 100°C for 0 min; 30°C /min to 200°C for 4 min, 15.083 min runtime |
| Inlet: | 220°C |
| Column Flow | 0.9mL/min |
| Gas: | Helium |
| Split: | 80:1 |
| Pressure: | 21.542 psi |
| Inlet: | Split/Split less |
| FID | 250°C, Hydrogen Flow 35.0mL/min, Air Flow 400.0mL/min, Constant Column and Makeup Flow 35mL/min |

Table 1: GC/FID Parameters

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

| Atomx Water Parameters | | | |
|--------------------------|-----------|---------------------------------|-----------------|
| Variable | Value | Variable | Value |
| Valve oven Temp | 140°C | Dry Purge Flow | 100mL/min |
| Transfer Line Temp | 140°C | Dry Purge Temp | 20°C |
| Sample Mount Temp | 90°C | Methanol Needle Rinse | Off |
| Water Heater Temp | 90°C | Methanol Needle Rinse Volume | 3.0mL |
| Sample Vial Temp | 20°C | Water Needle Rinse Volume | 7.0mL |
| Sample Equilibrate Time | 0.00 min | Sweep Needle Time | 0.50min |
| Soil Valve Temp | 50°C | Desorb Preheat Time | 245°C |
| Standby Flow | 10mL/min | GC Start Signal | Start of Desorb |
| Purge Ready Temp | 40°C | Desorb Time | 2.00 min |
| Condensate Ready Temp | 45°C | Drain Flow | 300mL/min |
| Presweep Time | 0.25 min | Desorb Temp | 250°C |
| Prime Sample Fill Volume | 3.0mL | Methanol Glass Rinse | Off |
| Sample Volume | 5.0mL | Number of Methanol Glass Rinses | 1 |
| Sweep Sample Time | 0.25 min | Methanol Glass Rinse Volume | 3.0mL |
| Sweep Sample Flow | 100mL/min | Number of Bake Rinses | 1 |
| Sparge Vessel Heater | On | Water Bake Rinse Volume | 7.0mL |
| Sparge Vessel Temp | 40°C | Bake Rinse Sweep Time | 0.25 min |
| Prepurge Time | 0.00 min | Bake Rinse Sweep Flow | 100mL/min |
| Prepurge Flow | 0mL/min | Bake Rinse Drain Time | 0.40 min |
| Purge Time | 11.00 min | Bake Time | 2.00 min |
| Purge Flow | 40mL/min | Bake Flow | 200mL/min |
| Purge Temp | 20°C | Bake Temp | 280°C |
| Condensate Purge Temp | 20°C | Condensate Bake Temp | 200°C |
| Dry Purge Time | 2.00 min | | |

Table 2: Atomx Water Parameters (Parameters highlighted in yellow were not used.)

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

| Atomx Soil Parameters | | | |
|-----------------------|-----------|------------------------------|------------------|
| Variable | Value | Variable | Value |
| Valve oven Temp | 140°C | Purge Time | 11.00 min |
| Transfer Line Temp | 140°C | Purge Flow | 40mL/min |
| Sample Mount Temp | 90°C | Purge Temp | 20°C |
| Water Heater Temp | 90°C | Condensate Purge Temp | 20°C |
| Sample Vial Temp | 40°C | Dry Purge Time | 2.00 min |
| Prepurge Time | 0.00 min | Dry Purge Flow | 100mL/min |
| Prepurge Flow | 0mL/min | Dry Purge Temp | 20°C |
| Preheat Mix Speed | Slow | Methanol Needle Rinse | Off |
| Sample Preheat Time | 0.00 min | Methanol Needle Rinse Volume | 3.0mL |
| Soil Valve Temp | 100°C | Water Needle Rinse Volume | 7.0mL |
| Standby Flow | 10mL/min | Sweep Needle Time | 0.25 min |
| Purge Ready Temp | 40°C | Desorbs Preheat Time | 245°C |
| Condensate Ready Temp | 45°C | GC Start Signal | Start of Desorbs |
| Presweep Time | 0.25 min | Desorbs Time | 2.00 min |
| Water Volume | 10mL | Drain Flow | 300mL/min |
| Sweep Water Time | 0.25 min | Desorbs Temp | 250°C |
| sweep Water Flow | 100mL/min | Bake Time | 2.00 min |
| Sparge Vessel Heater | Off | Bake Flow | 400mL/min |
| Sparge Vessel Temp | 20°C | Bake Temp | 270°C |
| Purge Mix Speed | Fast | Condensate Bake Temp | 200°C |

Table 3: Atomx Soil Parameters (Parameters highlighted in yellow were not used.)

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

| Atomx Methanol Extraction Parameters | | | |
|--------------------------------------|-----------|---------------------------------|----------------|
| Variable | Value | Variable | Value |
| Valve Oven Temp | 140°C | Dry Purge Flow | 100mL/min |
| Transfer Line Temp | 140°C | Dry Purge Temp | 20°C |
| Sample Mount Temp | 90°C | Methanol Needle Rinse | On |
| Soil Valve Temp | 100°C | Methanol Needle Rinse Volume | 2.0mL |
| Standby Flow | 10mL/min | Water Needle Rinse Volume | 7.0mL |
| Purge Ready Temp | 40°C | Sweep Needle Time | 0.25min |
| Condensate Ready Temp | 45°C | Desorb Preheat Temp | 245°C |
| Presweep Time | 0.25min | GC Start Signal | Star of Desorb |
| Methanol Volume | 7mL | Desorb Time | 2.00min |
| Sparge Vessel Heater | On | Drain Flow | 300mL/min |
| Sparge Vessel Temp | 40°C | Desorb Temp | 250°C |
| Prepurge Time | 0.00min | Methanol Glass Rinse | On |
| Prepurge Flow | 0mL/min | Number of Methanol Glass Rinses | 1 |
| Sample Mix Speed | Fast | Methanol Glass Rinse Volume | 3.0mL |
| Sample Mix Time | 4.00min | Number of Water Bake Rinses | 1 |
| Sample Settle Time | 2.00min | Water Bake Rinse Volume | 7.0mL |
| Sample Sweep Time | 0.25min | Bake Rinse Sweep Time | 0.25min |
| Sample Sweep Flow | 100mL/min | Bake Rinse Sweep Flow | 100mL/min |
| Purge Time | 11.00min | Bake Rinse Drain Time | 0.40min |
| Purge Flow | 40mL/min | Bake Time | 2.00min |
| Purge Temp | 20°C | Bake Flow | 200mL/min |
| Condensate Purge Temp | 20°C | Bake Temp | 280°C |
| Dry Purge Time | 2.00min | Condensate Bake Temp | 200°C |

Table 4: Atomx Methanol Extraction Parameters

Calibration / Results

A 50ppm working calibration standard was prepared in methanol using a 1000ppm Wisconsin GRO/PRVO standard. Calibration standards were then serially diluted with de-ionized water over a range of 5-200ppb. A 25ppm surrogate standard of fluorobenzene was prepared in methanol and transferred to one of the three standard addition vessels on the Atomx. The Atomx delivered the surrogate in 5µL aliquots to the sample to give a final concentration of 25ppb.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Using Agilent Chemstation software, the summation of the peak areas of each standard from methyl-tert-butyl-ether (MtBE) to naphthalene were used to generate the water and low-level soil calibration curves seen in **Figures 1-2**. The area of the surrogate fluorobenzene was subtracted from the total area to give the correct area required for the Wisconsin GRO calibration. The calibrations passed method criteria with a correlation coefficient (r^2) of 0.9978 for water and r^2 of 0.9990 for low level soil. Using the #11-VPH-Trap it will allow for the resolution of MtBE down to 5ppb allowing for a calibration curve to be run from 5 to 200ppb.




Figure 1: Wisconsin GRO water calibration on a range from 5-200ppb using #11-VPH-Trap on Teledyne Tekmar's Atomx




Figure 2: Wisconsin GRO low level soil calibration on a range from 5-200ppb using #11-VPH-Trap on Teledyne Tekmar's Atomx

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

The Wisconsin GRO method requires for Quality Control (QC) checks to demonstrate the capability to generate acceptable accuracy and precision.¹ For the water/soil, each matrix must run a RLCS, DLCS, LOD and LOQ. Each of these test listed above must fall between required limits that are set forth in the Wisconsin GRO method. Results of these test can be found in **Table 5-6**. **Table 5** shows the capability of the Atomx to generate acceptable accuracy and precision for both water and soil matrices using the #11-VPH-Trap. **Table 6** shows that the LOQs for both matrices are below the required limits set by the method.

| Water Quality Control | | | | |
|--|-------------------|--------------------|-----------|-------------------|
| | Theoretical (ppb) | Experimental (ppb) | %Recovery | Experimental %RSD |
| RLCS | 100 +/-20% | 95.73* | 95.73 | 2.47 |
| DLCS | 100+/-20% | 89.0 | 89.0 | |
| Low Level Soil Quality Control | | | | |
| RLCS | 100 +/-20% | 105.08* | 105.08 | 1.59 |
| DLCS | 100+/-20% | 97.48 | 97.48 | |
| High Soil Quality Control (Methanol Extracted) | | | | |
| | Theoretical (ppm) | Experimental (ppm) | %Recovery | Experimental %RSD |
| RLCS | 10+/-120-75% | 8.75* | 87.5 | 3.96 |
| DLCS | 10+/-120-75% | 8.61 | 86.1 | |

Table 5: Replicate Laboratory Control Spike (RLCS) and Duplicate Laboratory Control Spike (DLCS) Results for water and soil.
(*n= 5 replicates)

| Water/Soil Quality Control | | |
|----------------------------|------------|-----------|
| Matrix | LOD* (ppb) | LOQ*(ppb) |
| Water | 1.56 | 4.96 |
| Low Level Soil | 1.854 | 5.89 |
| High level Soil | 1.16 | 5.14 |

Table 6: Results for for the Lower Limit of Detection (LOD) and Limit of Quantitation LOQ for both water and soils.
(*n=7 replicates)

All the experimental values fall with in the acceptable range set forth by the Wisconsin GRO Method. The method states the LOQ for soil should be less then 10ppm and for groundwater it should be 0.1ppm or less.¹

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Velocity XPT (straight trap) with an AQUATek 100 Autosampler

Instrument Parameters

The second part of the #11-VPH-Trap study was to validate the trap on a system that requires a straight analytical trap. A Velocity XPT concentrator with an AQUATek 100 autosampler was used in tandem with a GC/FID. The Agilent 6890 GC/FID parameters can be found in **Table 1** while the Velocity XPT and AQUATek 100 instrument parameters can be found in **Table 7**.

| Velocity P&T and AQUATek 100 Parameters | | | |
|---|-----------------|-----------------------|------------|
| Variable | Value | Variable | Value |
| Valve Oven Temp | 150°C | Desorb Time | 2.00 min |
| Transfer Line Temp | 150°C | Desorb Temp | 250°C |
| Purge Ready Temp | 35°C | Drain Flow | 300 mL/min |
| Dry Flow Standby Temp | 40°C | Bake Time | 2.00 min |
| Standby Flow | 5 mL/min | Bake Temp | 280°C |
| Spurge Vessel Heater | On | Dry Flow Bake Temp | 270°C |
| Pre-Purge Time | 0.00 min | Bake Flow | 200 mL/min |
| Pre-Purge Flow | 40 mL/min | Pressurize Time | 0.35 min |
| Preheat Time | 1.00 min | Sample Transfer Time | 0.35 min |
| Sample Temp | 40°C | Rinse Loop Time | 0.30 min |
| Purge Temp | 20°C | Sweep Needle Time | 0.30 min |
| Purge Flow | 40 mL/min | Bake Rinse | On |
| Dry Purge Time | 2.00 min | Bake Rinse Cycle | 1 |
| Dry purge Temp | 20°C | Bake Rinse Drain Time | 0.35 min |
| Dry Purge Flow | 200 mL/min | Presweep Time | 0.25 min |
| GC Start | Start of Desorb | Water Temp | 90°C |
| Desorb Preheat Temp | 245°C | | |

Table 7: Velocity XPT and AQUATek 100 parameters. (AQUATek 100 parameters are listed in blue)

Calibration / Results

A 50ppm working calibration standard was prepared in methanol using a 1000ppm Wisconsin GRO/PRVO standard from Restek. Calibration standards were then serially diluted with de-ionized water over a range of 5-200ppb. A 25ppm surrogate of fluorobenzene was prepared in methanol and transferred to one of the two standard addition vessels on the AQUATek100. The autosampler delivered the surrogate in 5µL aliquots to the sample for a final concentration of 25ppb.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Using Agilent Chemstation software, the summations of each standard from methyl-tert-butyl-ether to naphthalene were used to generate the water calibration curves seen in **Figure 3**. The area of the surrogate fluorobenzene was subtracted from the total area to give the correct value required for the Wisconsin GRO calibration. The calibration pass method criteria with a correlation coefficient (r^2) of 0.9998.




Figure 3: Wisconsin GRO water calibration on a range from 5-200ppb using #11-VPH-Trap on Teledyne Tekmar's Velocity XPT Purge and Trap Concentrator AQUATEk 100 Autosampler.

Since the AQUATEk 100 is a water autosampler, only the water portion of the Wisconsin GRO method could be evaluated. RLCS, DLCS, LOD and, LOQ were performed. Results of these tests can be found in **Table 8-9**. **Table 7** shows the capability of the Velocity XPT with a AQUATEk 100 to generate acceptable accuracy and precision for water samples using the #11-VPH-Trap. **Table 9** shows that the LOQs for both matrices are below the required limits set by the method.

| Water Quality Control | | | | |
|-----------------------|-------------------|--------------------|-----------|-------------------|
| | Theoretical (ppb) | Experimental (ppb) | %Recovery | Experimental %RSD |
| RLCS | 100 +/-20% | 98.44* | 98.44 | 4.20 |
| DLCS | 100+/-20% | 97.10 | 97.10 | |

Table 8: Replicate Laboratory Control Spike (RLCS) and Duplicate Laboratory Control Spike (DLCS) Results for water. (*n= 5 replicates)

| Water Quality Control | | |
|-----------------------|------------|------------|
| Matrix | LOD* (ppb) | LOQ* (ppb) |
| Water | 1.22 | 3.90 |

Table 9: Results for for the LOD and LOQ for both water. (*n=7 replicates)

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

All the experimental values fall within the acceptable range set forth by the Wisconsin GRO Method. The method states for LOQ the groundwater should be 0.1ppm or less.¹

HT3™ Automated Headspace Vial Sampler

Instrument Parameters

The dynamic (trap) mode of the HT3™ headspace system was used to evaluate the GRO analysis. The HT3™ was connected to an Agilent 6890 GC/FID. HT3™ parameters are listed in **Table 10**. The Agilent 6890 GC/FID system parameters are listed in **Table 1**.

| Dynamic (Trap) | | | |
|---------------------|-----------|----------------|-----------|
| Variable | Value | Variable | Value |
| Platen/Sample Temp | 75°C | Dry Purge Time | 2.00 min |
| Valve Oven Temp | 150°C | Dry Purge Flow | 100mL/min |
| Transfer Line Temp | 150°C | Dry Purge Temp | 25°C |
| Standby Flow Rate | 50mL/min | Desorb Preheat | 245°C |
| Sample Preheat Time | 0.00 min | Desorb Temp | 250°C |
| Sweep Flow Rate | 200mL/min | Desorb Time | 1.00 min |
| Sweep Flow Time | 5.00 min | Trap Bake Temp | 270°C |
| Preheat Mixer | 0.00 min | Trap Bake Time | 3.00 min |
| Trap Material | #11 | Trap Bake Flow | 200mL/min |

Table 10: Dynamic HT3™ Headspace parameters.

Standard and Sample Preparation

A 100ppb stock Internal Standard (IS) solution was prepared by adding 50 µL of a 2000ppm fluorobenzene standard to 1 L of de-ionized water. Using a Wisconsin GRO/PVOC mix a 50ppm stock standard was made. Calibration standards were then serially diluted with de-ionized water containing a 100ppb IS of fluorobenzene. The working range of the calibration ranged from 5-200ppb for Wisconsin GRO. 5 mL of the working standards were transferred to 22 mL headspace vials with a glass syringe and capped for the GRO analysis.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Calibration / Results

The peak areas of the Wisconsin GRO standards from 5 to 200ppb were used to calculate the Response Factors (RF) by both external and internal standard calculations. The method requires %RSD less than 20% and correlation coefficients (r^2) greater than 0.99, these values are shown in **Table 11**.

| Calibration Data | | |
|------------------------|--------------------------------------|---|
| 5 mL Sample Data | %RSD of the Compound Response Factor | Correlation Coefficient of the Compound |
| Compound | #11-VPH-Trap | #11-VPH-Trap |
| MtBE | 8.59 | 0.9999 |
| Benzene | 3.92 | 0.9996 |
| Toluene | 2.32 | 0.9993 |
| Ethylbenzene | 2.50 | 0.9995 |
| m-, p-Xylene | 1.02 | 0.9995 |
| o-Xylene | 2.39 | 0.9998 |
| 1,3,5-Trimethylbenzene | 4.60 | 0.9999 |
| 1,2,4-Trimethylbenzene | 4.00 | 0.9999 |
| Naphthalene | 5.82 | 0.9998 |

Table 11: Calibration data showing the %RSD and correlation coefficient for the calibration curve (5-200ppb)

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent




Figure 4: Chromatogram showing the separation 100ppb Wisconsin GRO standard from Methanol (MeOH) on the Tekmar HT3™

Using the #11-VPH-Trap enabled resolution of MtBE from the methanol peak seen in **Figure 4**. For all ten compounds listed in **Table 11** all the %RSD were below 10%, while the correlation coefficient were all greater than the requirements of 0.99 or better.

Trap Comparison

Throughout this study, the #11-VPH-Trap was compared to other traps for differences in performing this method. The first difference is the ability of the #11-VPH-Trap to resolve MtBE down to 5ppb while the M-Trap can only resolve MtBE down to 10ppb. **Figure 5** shows an overlay comparing chromatograms of a 5ppb Wisconsin GRO water sample on the #11-VPH-Trap and the M-Trap.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent




Figure 5: 5ppb chromatogram showing an overlay of Wisconsin GRO standard on the #11-VPH-Trap (Blue) and M-Trap (Red).

Using the #11-VPH-Trap will allow the analyst to generate a calibration curve for all compounds over a range of 5-200ppb for the Atomx automated VOC sample prep system, the Velocity XPT concentrator, and the HT3™ automated headspace vial sampler.

The second difference uncovered in the trap comparison was each trap's methanol retention. If the trap retains too much methanol, the associated solvent peak could affect, or in some cases completely mask the MtBE due to the similar retention times. Figure 6 shows an overlay of the solvent peak of three traps: the #11-VPH-Trap, The #9 Trap, and the M-Trap.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent




Figure 6: An overlay of the solvent peaks for all three traps.

As **Figure 6** demonstrates the VOCARB 3000 Trap retains the most amount of methanol due to the fact that it is typically used in full VOC list which require stronger sorbent materials in order to trap lighter volatiles such as Freon for example. The #11-VPH-Trap and M-Trap are designed for the analysis of Wisconsin GRO and BTEX compounds, which do not contain many polar compounds. Since these two traps do not require full analyte list the sorbent beds can be modified so that each trap retains a smaller amounts of methanol. The methanol does not interfere with the first analyte of interest (MtBE), making the #11-VPH-Trap an idea analytical trap for GRO analysis.

The third and final difference while noticed between the #11-VPH-Trap and the M-Trap was the bake back pressure. The #11-VPH-Trap behaved similarly to the the Tekmar #9 trap, the M-Trap showed evidence of increasing bake back pressure. This increasing pressure can lead to early degradation of the trap resulting in more frequent replacement of the analytical trap.

Evaluation of a New Analytical Trap for Gasoline Range Organics Analysis

Application Note

By Tyler Trent

Conclusion

The #11-VPH-Trap is the latest advancement in analytical traps by Teledyne Tekmar. This analytical trap provides the features of the Tekmar #9 analytical trap, but also reduces the sorbent bed so that lesser amount of methanol is retained. The trap keeps its ability to trap the polar and non-polar compounds found in gasoline range organics (ie: Wisconsin GRO and BTEX methods).

This study validates the use of the #11-VPH-Trap in three of Tekmar products the Atomx, Velocity XPT with an AQUATek 100 autosampler, and the HT3™ headspace analyzer. By using the #11-VPH-Trap the solvent peak is greatly reduced and does not interfere with the first analyte, MtBE. The #11-VPH-Trap also has a reduced bake pressure than the M-Trap which translates to a longer life cycle of the trap and reduces overall operation costs.

References

1. Modified GRO Method for Determine Gasoline Range Organics Wisconsin DNR September 1995 PUBL-SW-140